



Certificate of Analysis

Cat. No.	HMD002101
Enzyme Name	T4 Gene 32 Protein
Batch Number	HH20220126
Storage Temperature	-20°C
Recommended re-test date	2023.01.25
Sources	Recombinant microorganism

Test items	Specifications	Results
Endonuclease Activity	Not detectable	PASS
Exonuclease Activity	Not detectable	PASS
Nicking Activity	Not detectable	PASS
RNase Activity	Not detectable	PASS
<i>E.coli</i> DNA	≤1 copy/10 μg	PASS
Protein Concentration (A280)	10 mg/ml ± 5%	PASS
Purity	≥95%	PASS

Analyzed by

Authorized by

Report Date

